Table of Contents

ENGLISH IFU.....	2
SPANISH IFU.....	6
FRENCH IFU	10
PORTUGUESE (BRAZIL) IFU	14
GERMAN IFU.....	16
DANISH IFU.....	22
GREEK IFU.....	26
CROATIAN IFU	30
ITALIAN IFU.....	34
DUTCH IFU.....	38
NORWEGIAN IFU.....	42
POLISH IFU	46
PORTUGUESE IFU	50
ROMANIAN IFU	54
SWEDISH IFU.....	58
SLOVENIAN IFU	62
TURKISH IFU.....	66

Instructions for Use

COPYRIGHT© 2024 BY SEGAMI CORPORATION, ALL RIGHTS RESERVED. Processing software and the related documentation contain confidential and/or proprietary information which is protected by copyright and international treaty provisions. Any unauthorized use, disclosure, assignment, transfer or reproduction of this confidential information will be prosecuted to the full extent of the Law.

Published by Segami Corporation, 8310 Guilford Rd, Suite A, Columbia, MD USA 21046-2818

Regulatory Standards

The following standards have been applied in development of this product: ISO 13485:2016/EN ISO 13485:2016 International Standard, Medical Device Directive (MDD) 93/42/EEC EU Medical Device Regulations 2017/745 in transition, Health Canada Medical Device Regulations (CMDR) SOR/98-282, ANVISA (Brazil) Technical Regulation for Good Manufacturing Practices of Medical Devices, and the FDA's Quality System Regulations (QSRs). The last three countries are three member states (Canada, Brazil, and the US) of the Medical Device Single Audit Program (MDSAP).

CE Conformity

This product is classified as Class IIa. Segami Corporation declares that this product meets the provision of the Council Directive 2007/47/EC for Medical Devices, and therefore bears the CE mark of conformity (Figure A). For the marking of Segami's Authorized Representative, see Figure B. For Segami's Importer symbol, see Figure C. These markings also appear in the Oasis application, under Help - About.

This product is manufactured by:

Segami Corporation

8310 Guilford Road

Suite A

Columbia, MD, USA 21046

+1-410-381-2311 – phone

+1-410-381-2320 – fax

Figure A

The authorized CE Representative for Segami Corporation is:

MedEnvoy Global B.V.

Prinses Margrietplantsoen 33 – Suite 123
2595 AM The Hague
The Netherlands

Figure B

The Importer (Economic Operator) for EU is:

MedEnvoy Global B.V.

Prinses Margrietplantsoen 33 – Suite 123
2595 AM The Hague
The Netherlands


Figure C

Program Description

Oasis is an advanced visualization, quantification, and reporting software solution for nuclear medicine diagnostic imaging. It integrates multiple modalities — SPECT, PET, and CT — for the accurate interpretation of organ function and anatomy in a single, comprehensive environment. The intuitive, user-centered environment guides the user while providing configurable preferences for workflow personalization. The core functionality is to provide an accurate environment for quantification and display of nuclear medicine images and related other modalities. It supports a wide variety of study types, based on European and North American published standards and guidelines. Oasis provides a non-manufacturer specific solution for review of images from a wide variety of available scintigraphic cameras. Data may be passed to Oasis either from a connected or integrated Oasis server/database or from various commercially available PACS system databases. The system is intended for use in imaging centers and hospitals by qualified nuclear medicine personnel trained in the use of the Oasis product.

Users without normal color vision should use select non-quantitative greyscale color tables and may require additional training.

Indications for Use

Oasis is indicated for processing and review of diagnostic medical images produced by DICOM based multi-modality sources; in other words, it provides the physicians with a tool to get, format, store and display diagnostic image data in all traditional manners; to add three-dimensional and multi-planar renderings to the displays; to fuse multimodality images and to display certain kinetic attributes of the imaged organs and systems.

Analysis in Oasis is NOT intended to be used as a replacement for visual interpretation nor as a diagnostic tool without other clinical and laboratory information.

Oasis is not used to treat the patient, to monitor vital signs or to give a direct diagnosis of disease.

Intended User Profile

The intended Oasis users are trained in the use of this product and are qualified nuclear medicine or qualified Radiologists. For users not trained in the use of Oasis, please contact support@segamicorp.com.

Residual Risks

Segami Corporation complies with EN ISO 14971:2019 and has both safety and security risk management programs. Known and foreseeable clinical risks have been considered for Oasis in accordance with Doc#:563 QPR-4.1 Risk Management and Doc#:709 QPR 4.3 Cybersecurity Policy. Segami has determined there are no residual risks or undesirable side-effects (e.g., adverse effects, complications) to be conveyed to the patient or user via the IFU, though the importance of being properly trained in the use of Oasis is stressed. Users not adequately trained should not use the product and could contribute to the misuse or off-label use of Oasis which can lead to increased risk for both safety and security.

Security Recommendations

Since the Oasis client/server software is installed on customer provided hardware within the site's network, it is expected that these institutions and medical clinics have established cybersecurity policies to protect against threats to their technical infrastructure. To request Doc#:434 Hardware Specifications, please contact support@segamicorp.com.

Effective Information Assurance/Security policies should consider the following:

People – Employees are the heart of information systems. Administrators and users of information systems must follow policies and practice for designing a good system. They must be informed regularly regarding information system and ready to act appropriately to safeguard system.

Policy & Practice – Organizations should define rules in the form of policies that every employee must follow. These policies must be practiced in order to properly handle sensitive information whenever system gets compromised.

Technology – Appropriate technology such as firewalls, routers, and intrusion detection must be used in order to defend systems from vulnerabilities and threats. The technology used must facilitate quick response whenever information security gets compromised.

Using Oasis

Opening Oasis

- Double-click the Oasis icon on the desktop.
- Once the login screen appears, enter the username and password and click [Login]. The patient/study browser screen will appear.

Closing Oasis

- First, close any application currently in use by clicking the [Return to Browser] button.
- Then close the patient/study browser by selecting the File ► Logout menu or by clicking 'X' in the top right corner.

About Oasis Display

- Because the Oasis display is customized automatically for each protocol, the number of screens, pages, and viewports rendered will vary, depending upon the clinical application protocol in use.

Access to Oasis Server Web Services

- When successfully logged into Oasis client select Tools ► Administrative Tools to access the Oasis Server interface to perform operations such as Query/Retrieve, image or study deletion, and other data administration tasks.

Location of Help Files

- Press [F1] to display a list of keyboard shortcuts. Press  to view the User Manual.

Support

The product is distributed and supported only by Segami or its qualified sales representatives. Please contact your local sales representative. Email support@segamicorp.com to request support. To provide remote access to Segami, select Help ► Remote Support to activate direct support from Segami qualified personnel.

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COPYRIGHT ©2024 PAR SEGAMI CORPORATION, TOUS DROITS RÉSERVÉS
COPYRIGHT© 2024 SEGAMI CORPORATION, ALLE RECHTE VORBEHALTEN. Das Computerprogramm und die zugehörige Dokumentation enthalten vertrauliche und/oder geschützte Informationen, die durch Urheberrechte und internationale
COPYRIGHT© 2024 SEGAMI CORPORATION, SVA PRAVA PRIDRŽANA. Softver za obradu i pripadajuća dokumentacija sadrže povjerljive i/ili zaštićene podatke koji su zaštićeni autorskim pravima i odredbama međunarodnih ugovora. Svako neovlašteno korištenje, otkrivanje, ustupanje, prijenos ili reprodukcija ovih povjerljivih informacija kazneno će se goniti u punoj mjeri Zakona.
COPYRIGHT © 2024 AV SEGAMI CORPORATION, ALLE RETTIGHETER FORBEHOLDES. Programvaren og tilhørende dokumentasjon inneholder konfidensiell og/eller proprietær informasjon som er beskyttet av opphavsrett og internasjonale
COPYRIGHT© 2024 BY SEGAMI CORPORATION, WSZELKIE PRAWA ZASTRZEŻONE. Oprogramowanie przetwarzające i związana z nim dokumentacja zawierają informacje poufne i/lub
COPYRIGHT© 2024 DA SEGAMI CORPORATION, TODOS OS DIREITOS RESERVADOS. O software de processamento e a documentação relacionada contém informações confidenciais e/ou proprietárias protegidas por direitos autorais